

clinell®

www.clinell.com

ALCOHOL WIPES

Surface disinfection

Alcohol wipes are suitable for surfaces of non-invasive medical devices.

Kills 99.999% of germs

Includes MRSA, *E. coli*,
Pseudomonas, *Enterococcus*,
Candida and *Aspergillus*
according to EN1276, EN13727,
EN14561, EN13624 and EN14562.

The pharmaceutical grade isopropyl alcohol evaporates and is suitable for surfaces such as glass, stainless steel and monitor screens.



ALCOHOL WIPES

Disinfection wipes for surface disinfection of non-invasive medical devices.



Alcohol Wipes

Clinell Alcohol Wipes are high quality durable wipes saturated with pharmaceutical grade 70% isopropyl alcohol. Our wipes act as a rapid disinfectant for surfaces and equipment, with a proven bactericidal action.

Our flow wrap wipes contain 25gsm nonwoven spun bond.

Unique packaging

Brings alcohol wipe packaging into the 21st century - wipes are dispensed one at a time, so the top of the next wipe is not dry which commonly occurs with tub dispensers. This eliminates the infection control risk of dry contaminated wipes. They are also environmentally friendly as there is much less packaging to dispose of compared to standard alcohol tubs.

FOR USE ON THE FOLLOWING SURFACES AND EQUIPMENT

- Door handles
- Bed frames
- Trolleys
- Walking frames
- Wheelchairs
- Radiology equipment
- Monitors
- Blood pressure machines & cuffs
- Thermometers
- Keyboards & accessories
- Telephones
- Stethoscopes
- ECG machines
- Ultrasound probes
- Doppler probes
- Pulse oximeter probes

BACTERIA	TEST
<i>Escherichia coli (E. coli)</i>	EN 1276
<i>Enterococcus hirae</i>	EN 1276
<i>Enterococcus hirae</i>	EN 14561
<i>Enterococcus hirae</i>	EN 13727
<i>Staphylococcus aureus</i>	EN 1276
<i>Staphylococcus aureus</i>	EN 14561
<i>Staphylococcus aureus</i>	EN 13727
<i>Pseudomonas aeruginosa</i>	EN 1276
<i>Pseudomonas aeruginosa</i>	EN 14561
<i>Pseudomonas aeruginosa</i>	EN 13727
FUNGI	TEST
<i>Candida albicans</i>	EN 14562
<i>Candida albicans</i>	EN 13624
<i>Aspergillus niger</i>	EN 14562
<i>Aspergillus niger</i>	EN 13624

Kills at least
99.999% of germs.

PRODUCT	UNIT OF ISSUE	ORDER CODE	NHSSC
Alcohol Wipes Small 100	Pack of 100	CAW100S	VJT428
Alcohol Wipes Large 100	Pack of 100	CAW100L	VJT421
Alcohol Wipes Large 150	Pack of 150	CAW150L	VJT448

USE DISINFECTANTS SAFELY. ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE. ALWAYS FOLLOW THE MEDICAL DEVICE MANUFACTURER'S CLEANING AND DISINFECTION GUIDELINES PRIOR TO USE WITH CLINELL ALCOHOL WIPES.

CE 0050 CLASS IIa MEDICAL DEVICE

Conforms to the Medical Device Directive, 93/42/EEC.

JBN200318